What is the N-back Test?

The **N-back Test** is a **working memory task** widely used in cognitive neuroscience to measure **mental workload (MW)**. It requires participants to **remember a sequence of items** and respond if the current item matches the one presented **N steps earlier**.

Purpose in This Dataset

To induce and measure three distinct levels of mental workload:

- 1-back (Low Load): Easy memory recall.
- 1-back Arithmetic (Medium Load): Basic math + memory.
- 2-back Dual Task (High Load): Combines memory + spatial + math tasks.

How It Was Used in the Experiment

Participants:

- 16 male participants
- Aged 20 to 60
- Scientists or researchers from universities (students, PhDs, professors)

Session Timeline:

Each subject completed:

- Baseline phase (10 min): Watched relaxing video (no task)
- Task phase (20 min): Played the N-back game
- Recovery phase (10 min): Resting state

Each subject performed all 3 N-back test variants. The order was randomized to avoid sequence bias.

Data Files and Structure

Located in: data_n_back_test/

1. eeg.parquet

Contains the raw EEG signals, plus:

- 139 columns, including:
 - datetime, timestamp
 - o **14 EEG channel values** (e.g., AF3, AF4, T7, T8, etc.)
 - o **5 power band values** (theta, alpha, beta low, beta high, gamma)
 - Quality metrics (signal quality scores per channel)
- Plus 3 metadata columns:

```
subject: e.g., subject_01
```

test: 1 (low), 2 (medium), 3 (high)

phase: 1 (baseline), 2 (task), 3 (recovery)

2. game_scores.parquet

Contains:

- Game performance scores for each subject and test
- Since each test runs for 20 minutes with multiple trials, scores are stored as a list of correct responses per trial
- Metadata:

- subject
- o test

3. tlx_answers.parquet

Contains:

- NASA-TLX questionnaire results for each test and subject
 - Measures perceived mental demand, effort, frustration, etc.
- Metadata:
 - subject
 - o test

Types of N-back Tasks Used

- 1-Back Position Test (Low Workload)
 - A square appears in one of 9 grid positions.
 - The subject presses a key if the current square **matches the previous one**.
 - Emphasizes visual-spatial memory.
- 1-Back Arithmetic Test (Medium Workload)
 - A number appears on the screen (0–9).
 - An arithmetic operation (e.g., +, -, ×, ÷) is heard.
 - The subject computes: current op previous and types the result.

2-Back Dual Task (High Workload)

- Number appears in a position on a 3×3 grid.
- The subject must:
 - Press a key if the **position matches** that from 2 steps before.
 - Solve an **arithmetic operation** using the number from 2 steps ago.

This test imposes high **cognitive load**, combining **working memory**, **attention**, **math**, **and spatial reasoning**.

```
['timestamp', 'EEG.Counter', 'EEG.Interpolated', 'EEG.AF3', 'EEG.F7',
'EEG.F3', 'EEG.FC5', 'EEG.T7', 'EEG.P7', 'EEG.01', 'EEG.02', 'EEG.P8',
'EEG.T8', 'EEG.FC6', 'EEG.F4', 'EEG.F8', 'EEG.AF4', 'EEG.RawCq',
'EEG.Battery', 'MarkerIndex', 'MarkerType', 'MarkerValueInt',
'EEG.MarkerHardware', 'CQ.AF3', 'CQ.F7', 'CQ.F3', 'CQ.FC5', 'CQ.T7',
'CQ.P7', 'CQ.01', 'CQ.02', 'CQ.P8', 'CQ.T8', 'CQ.FC6', 'CQ.F4', 'CQ.F8',
'CQ.AF4', 'PM.Engagement.IsActive', 'PM.Engagement.Scaled',
'PM.Engagement.Raw', 'PM.Engagement.Min', 'PM.Engagement.Max',
'PM.Excitement.IsActive', 'PM.Excitement.Scaled', 'PM.Excitement.Raw',
'PM.Excitement.Min', 'PM.Excitement.Max', 'PM.LongTermExcitement',
'PM.Stress.IsActive', 'PM.Stress.Scaled', 'PM.Stress.Raw',
'PM.Stress.Min', 'PM.Stress.Max', 'PM.Relaxation.IsActive',
'PM.Relaxation.Scaled', 'PM.Relaxation.Raw', 'PM.Relaxation.Min',
'PM.Relaxation.Max', 'PM.Interest.IsActive', 'PM.Interest.Scaled',
'PM.Interest.Raw', 'PM.Interest.Min', 'PM.Interest.Max',
'PM.Focus.IsActive', 'PM.Focus.Scaled', 'PM.Focus.Raw', 'PM.Focus.Min',
'PM.Focus.Max', 'POW.AF3.Theta', 'POW.AF3.Alpha', 'POW.AF3.BetaL',
'POW.AF3.BetaH', 'POW.AF3.Gamma', 'POW.F7.Theta', 'POW.F7.Alpha',
'POW.F7.BetaL', 'POW.F7.BetaH', 'POW.F7.Gamma', 'POW.F3.Theta',
'POW.F3.Alpha', 'POW.F3.BetaL', 'POW.F3.BetaH', 'POW.F3.Gamma',
'POW.FC5.Theta', 'POW.FC5.Alpha', 'POW.FC5.BetaL', 'POW.FC5.BetaH',
'POW.FC5.Gamma', 'POW.T7.Theta', 'POW.T7.Alpha', 'POW.T7.BetaL',
'POW.T7.BetaH', 'POW.T7.Gamma', 'POW.P7.Theta', 'POW.P7.Alpha',
'POW.P7.BetaL', 'POW.P7.BetaH', 'POW.P7.Gamma', 'POW.O1.Theta',
'POW.O1.Alpha', 'POW.O1.BetaL', 'POW.O1.BetaH', 'POW.O1.Gamma',
'POW.02.Theta', 'POW.02.Alpha', 'POW.02.BetaL', 'POW.02.BetaH',
```

```
'POW.O2.Gamma', 'POW.P8.Theta', 'POW.P8.Alpha', 'POW.P8.BetaL',
'POW.P8.BetaH', 'POW.P8.Gamma', 'POW.T8.Theta', 'POW.T8.Alpha',
'POW.T8.BetaL', 'POW.T8.BetaH', 'POW.T8.Gamma', 'POW.FC6.Theta',
'POW.FC6.Alpha', 'POW.FC6.BetaL', 'POW.FC6.BetaH', 'POW.FC6.Gamma',
'POW.F4.Theta', 'POW.F4.Alpha', 'POW.F4.BetaL', 'POW.F4.BetaH',
'POW.F4.Gamma', 'POW.F8.Theta', 'POW.F8.Alpha', 'POW.F8.BetaL',
'POW.F8.BetaH', 'POW.F8.Gamma', 'POW.AF4.Theta', 'POW.AF4.Alpha',
'POW.AF4.BetaL', 'POW.AF4.BetaH', 'POW.AF4.Gamma', 'datetime', 'subject',
'test', 'phase']
```

EEG File Column Breakdown

1. Time & Metadata

Column	Description
timesta mp	Unix time or elapsed time (in seconds or ms).
datetim e	Actual date & time of recording.
subject	Subject ID, e.g., subject_01.
test	N-back test type: 1 = Low, 2 = Medium, 3 = High workload.
phase	1 = Baseline, 2 = Task, 3 = Recovery.

2. Raw EEG Signals

14 EEG electrodes as per the **10-20 system**:

Column	Electrode Location	
EEG.AF3, EEG.AF4	Anterior frontal	
EEG.F3, EEG.F4, EEG.F7, EEG.F8	Frontal	
EEG.FC5, EEG.FC6	Fronto-central	

EEG.T7, EEG.T8	Temporal
EEG.P7, EEG.P8	Parietal
EEG.01, EEG.02	Occipital

Other columns:

• EEG.Counter, EEG.Interpolated, EEG.RawCq, EEG.Battery → device-specific telemetry data.

3. Event Markers (Optional)

Column	Description
--------	-------------

MarkerIndex, MarkerType, MarkerValueInt, EEG.MarkerHardware

Useful if task events (stimuli) were annotated.

4. Channel Quality (CQ)

These represent signal quality per electrode.

Example	Description
---------	-------------

CQ.AF3, CQ.F3, ..., Channel quality scores (typically 0–4 or a percentage).

These help you detect **bad or noisy channels** during preprocessing.

5. Performance Metrics (PM)

Category Examples

EngagementPM.Engagement.Raw, .Scaled, .IsActive, etc.

Excitement, Stress, Relaxation, Interest, Focus

Similar structure: each with .IsActive, .Raw, .Scaled, .Min, .Max.

These are **Emotiv's proprietary real-time affective state indicators**, inferred from EEG. They **aren't perfect**, but they can offer high-level estimates of cognitive states like **focus and stress**.

6. Power Spectral Features (Frequency Bands)

These are the **core EEG workload features** you'll want to use for ML or statistical analysis.

For **each EEG channel**, power values are computed for:

- Theta (4–8 Hz)
- Alpha (8–12 Hz)
- BetaL (12–18 Hz)
- BetaH (18-25 Hz)
- Gamma (>25 Hz)

For example:

- POW.AF3.Theta → Theta power at AF3.
- POW.F7.BetaH → Beta-high power at F7.

There are 14 electrodes × 5 frequency bands = 70 spectral features.

14 EEG channels:

- Frontal: AF3, AF4, F3, F4, F7, F8
- Central: FC5, FC6
- Temporal: T7, T8

- Parietal: P7, P8
- Occipital: 01, 02

Metadata:

- datetime: Timestamp of each sample
- subject: Always subject_01 in this case
- test: N-back condition \rightarrow 1 = low, 2 = medium, 3 = high
- phase: 1 = baseline, 2 = task, 3 = recovery

Mental Workload EEG Dataset – N-back Test Focus

This dataset is part of the broader research:

"EEG Dataset Collection for Mental Workload Predictions in Flight-Deck Environment"

GENERAL DATASET OVERVIEW

Recording Devices EEG: Emotiv Epoc X (14 channels, 128 Hz), ECG (Shimmer3,

Suunto Ambit3)

Experiments Conducted

3: N-back Test 🧠, Heat-the-Chair Game 🎮, Flight Simulator 🛫

Subjects 16 for N-back, 17 for Chair game, 2 pilots for flights

Data Type Raw EEG signals + Power Band data + Game scores + Subjective

TLX ratings

Sampling Rate 128 Hz (EEG), various physiological metrics per second (ECG,

power bands)

Data Format .parquet, .csv, .json

Labels Theoretical workload level (test), self-perceived TLX rating

Use License CC BY-SA

Link Dataset DOI

FOCUS: N-back Test (Memory Load Assessment)

Objective

To study brain activity and mental workload under controlled memory tasks of increasing difficulty using EEG.

Participants

- **16 male subjects** (aged 20–60)
- From 3 university labs (students, PhD researchers, professors)

N-back Test Structure

Each participant completed 3 test variants (randomized order):

Tes t	Description	Workload Level
1	1-back Position: Match square pos.	Low
2	1-back Arithmetic: Solve math op.	Medium
3	2-back Dual: Math + Position (2-back)	High

Session Timeline Per Test

Phase	Duration	Description
Baseline	10 min	Subject watches calming video
Task	20 min	Performs the N-back task
Recovery	10 min	Rest phase after task
Total	40 min	\rightarrow ~307,200 samples/test @ 128 Hz

Each subject = 3 test sessions → ~921,600 expected samples

III EEG Data Details

EEG Channels

14 channels based on the 10-20 system:

AF3, AF4, F3, F4, F7, F8, FC5, FC6, T7, T8, P7, P8, O1, O2

Recorded Data Includes:

- Raw EEG voltages (in μV)
- EEG channel quality
- Power spectral features:
 - Theta (4–8 Hz)
 - Alpha (8–12 Hz)
 - Beta Low (12–18 Hz)
 - Beta High (18–25 Hz)
 - o Gamma (>25 Hz)

- Mental state metrics: focus, stress, engagement, etc.
- Timestamps and phases
- Game performance scores
- NASA-TLX guestionnaire scores (self-perceived workload)

Workload Labels & Validation

Each test is labeled:

- test = 1 → low workload
- test = 2 → medium workload
- test = 3 → high workload

Validation was done using:

- 1. Game performance (decreases with load)
- 2. TLX scores (increase with load)
- 3. EEG power peaks (band-wise differences by test/phase)
- 4. Deep learning classifiers (CNNs trained on data)

Files Involved (N-back Test Folder)

bash

CopyEdit

```
\begin{array}{lll} data\_n\_back\_test/ \\ & \longmapsto & eeg.parquet & \to & Raw \ EEG \ + \ power \ + \ metadata \\ & \longmapsto & game\_scores.parquet & \to & Scores \ for \ each \ subject/test \\ & \longmapsto & tlx\_answers.parquet & \to & Subjective \ workload \ (NASA-TLX) \end{array}
```

```
Full Dataset Shape (rows, columns): (15294488, 18)
(i) Global Start Time: 2020-01-15 08:10:42.007477
Global End Time: 2020-11-19 13:48:46.993748
Total Duration: 309 days 05:38:04.986271
Total Unique Subjects: 16
Subjects: ['subject_01' 'subject_02' 'subject_03' 'subject_04' 'subject_05'
'subject 06' 'subject 07' 'subject 08' 'subject 09' 'subject 10'
'subject_11' 'subject_12' 'subject_13' 'subject_14' 'subject_15'
'subject 16']
Row count per subject:
subject
subject 01
            974072
subject_02
            973157
subject 03
            959145
subject_04
            998694
subject 05
            983382
subject_06
            978269
subject 07
            950255
subject 08
            994502
subject 09
            878821
subject 10
            836755
subject_11
            1047537
subject_12
            983090
subject 13
            976636
subject 14
            964373
subject_15
            831361
subject_16
            964439
Name: count, dtype: int64
Rows per phase:
```

phase

- 1 3964087
- 2 7694147
- 3 3636254

Name: count, dtype: int64

Rows per test:

test

- 1 5006127
- 2 5095468
- 3 5192893

Name: count, dtype: int64

```
Rows per (subject, test, phase):
    subject test phase row_count
0 subject_01
               1
                   1
                       82609
                   2
1 subject 01
                       168144
               1
2 subject 01
                       80959
3 subject_01
               2
                   1
                       80581
4 subject_01
               2
                   2
                       159001
138 subject_16 2
                    2
                       158240
139 subject 16
                2
                    3 80713
140 subject_16
                3
                    1 81853
                    2 165093
141 subject 16
                3
142 subject_16
                        80834
```

[143 rows x 4 columns]

Subject 1

- Data shape (rows, columns): (974072, 18)

 Start time: 2020-01-17 09:57:02.001394

 End time: 2020-01-17 17:17:48.994743
- Total duration: 0 days 07:20:46.993349
- Rows per phase:

phase

- 1 244917
- 2 487420
- 3 241735

Name: count, dtype: int64

Rows per test:

test

- 1 331712
- 2 320163
- 3 322197

Name: count, dtype: int64

- Rows per (test, phase): test phase row_count
- 0 1 1 82609
- 1 1 2 168144
- 2 1 3 80959
- 3 2 1 80581
- 4 2 2 159001
- 5 2 3 80581

6 3 1 81727 7 3 2 160275

8 3 3 80195